



***LIGO Laboratory / LIGO Scientific Collaboration***

T1900395-v1

*LIGO*

June 2019 5

**Incoming Test of SQZ Spare Mephisto 1000NE Laser**

Jason Oberling

Distribution of this document:

LIGO Scientific Collaboration

This is an internal working note  
of the LIGO Laboratory.

**California Institute of Technology**  
**LIGO Project, MS 100-36**  
**1200 E. California Blvd.**  
**Pasadena, CA 91125**  
Phone: (626) 395-2129  
FAX: (626) 304-9834  
E-mail: [info@ligo.caltech.edu](mailto:info@ligo.caltech.edu)

**Massachusetts Institute of Technology**  
**LIGO Project, NW22-295**  
**185 Albany St.**  
**Cambridge, MA 02139**  
Phone: (617) 253-4824  
FAX: (617) 253-7014  
E-mail: [info@ligo.mit.edu](mailto:info@ligo.mit.edu)

**LIGO Hanford Observatory**  
**P.O. Box 159**  
**Richland, WA 99352**  
Phone: (509) 372-8106  
FAX: (509) 372-8137

**LIGO Livingston Observatory**  
**P.O. Box 940**  
**Livingston, LA 70754**  
Phone: (225) 686-3100  
FAX: (225) 686-7189

## 1. Scope and Purpose

This document shows the results of the incoming tests performed on the Coherent Mephisto 1000NE S/N 5834 laser system, hereafter referred to as the NPRO. This is the spare NPRO laser for the SQZ subsystem.

## 2. Introduction

Upon unpacking, the NPRO powered up without issue. Below are the serial numbers of the NPRO and its power supply, as well as a list of equipment used for the testing.

### 2.1 Serial Numbers

- Mephisto 1000NE, P/N 1309412, S/N GDP.1309412.5834
- Mephisto Controller, P/N 1309412, S/N GDP.1309412.5834

### 2.2 Equipment

- Power meter: Ophir Vega, P/N 7Z01560, S/N 730491
- Calorimeter: Ophir 10A-V2-SH, P/N 1Z02146, S/N 75126
- Beam profiler: Thorlabs BP209-VIS
- Photodiode: Thorlabs PDA36A
- Signal Analyzer: SRS SR785

## 3. Laser Settings

The data sheet for the NPRO is located on the DCC file card for this test report. For completeness, the NPRO settings are shown in Table 1.

**Table 1: NPRO Settings**


NPRO crystal temperature	26.20 °C
Diode current	2.14 A
Diode A temperature	26.01 °C

## 4. Output Power

The output power was measured to be 1.12 W after a warm-up period of ~5 minutes.

## 5. Free-space Beam Propagation

The free-space beam propagation was measured with the beam profiler listed above and is shown in Figure 1 along with a Gaussian fit. The transmitted beam through a 98/2 (%R/%T) beamsplitter was used to do this measurement. For the purpose of this measurement, zero is the flat front face of the NPRO laser housing.



**Figure 1:** The measured free-space beam propagation.

Table 2 summarizes the NPRO beam waist radius and location. As stated above, zero is the flat front face of the NPRO laser head housing.

**Table 2:** NPRO Waist Radius and Location

	Waist Radius ( $\mu\text{m}$ )	Waist Location (mm)
<b>Horizontal</b>	146	-89.4
<b>Vertical</b>	157	-94.4

## 6. Relative Power Noise

The measured low-frequency relative intensity noise (RIN) is shown in Figure 2; the measurement was performed with a SR785 Digital Signal Analyzer from Stanford Research Systems. It should be noted that Coherent does not give a performance specification for the RIN below 10 kHz. It simply lists the RIN as better than -140 dB/Hz at frequencies greater than 10 kHz, which the NPRO satisfies.




Figure 2: Measured RIN of Mephisto 1000NE SN5834